

**Stamatatos Lab**  
**ELISA Protocol**

Contact: Zane Kraft at [zane.kraft@seattlebiomed.org](mailto:zane.kraft@seattlebiomed.org)

1. 50ng of protein in 100  $\mu$ L of NaHCO<sub>3</sub> (100mM, pH 9.0) is added to each well of a high binding 96 well plate, incubate O.N. at RT
2. Block the plate with 250ul/well of 3% BSA in PBS, incubate 2hrs at RT
3. Serial dilute the Mab or Sera in PBS + 1% BSA + 0.03% Tween20 (dilution buffer) on the blocked ELISA plate, incubate 1hr at 37°C
4. Wash the plate 4 times
5. Dilute the HRP conjugated secondary Ab 1:3,000 in dilution buffer and add 100ul to each well, incubate 1hr at 37°C
6. Repeat the wash step
7. Add 50 $\mu$ L per well of TMB substrate, incubate 3 min at RT
8. Stop development by adding an equal volume of 1N H<sub>2</sub>SO<sub>4</sub>
9. Absorbance is measured at 450 nm

To determine the avidity: sera diluted to 1:500 is added to the blocked plate at step 3 and serially diluted ammonium thiocyanate, ranging from 10M to 1M, is then added to the diluted sera and incubated at RT for 20 min on a plate shaker