



# **INSTRUCTION MANUAL**

---

## **FAS FACIAL PHOTOGRAPHIC ANALYSIS SOFTWARE**

**VERSION 1.0.0**

**FAS DIAGNOSTIC AND PREVENTION NETWORK  
FAS DPN  
UNIVERSITY OF WASHINGTON  
SEATTLE WASHINGTON**

**SUSAN ASTLEY, PH.D.  
2003**

## Table of Contents

	Acknowledgments.....	iii
	License Agreement .....	iv
<b>PART I.</b>	<b>Introduction</b>	
	What is FAS.....	1
	The FAS 4-Digit Diagnostic Code.....	1
	Facial Features of FAS.....	2
	FAS Facial Analysis Software.....	2
<b>PART II.</b>	<b>Getting Started</b>	
	Professional Qualifications for Use of this Software.....	3
	Computer System.....	3
	Minimum Requirements .....	3
	Recommended Requirements .....	3
	To Install the Facial Software .....	3
	Help.....	3
	To Open the Software .....	4
	To Open the “John Doe” Demonstration Case .....	4
<b>PART III</b>	<b>Taking the Facial Photographs</b>	
	How to Take the 3 Standardized Facial Photographs .....	5
	Internal Measure of Scale .....	5
	Frontal, ¾, Lateral Pictures.....	5
	Facial Expression.....	5
	Rotation.....	6
	Additional Training / Instructional Resources.....	6
	Camera / Image Resolution Specifications	
	Digital Camera.....	7
	Film Camera.....	7
	Image Resolution .....	7
	Image Formats .....	7
<b>PART IV.</b>	<b>USING THE MEASUREMENT TOOLS, PICTORIAL GUIDES AND SCORING TABLES</b>	
	Tools and Guides to Measure Facial Features	
	Tools .....	8
	Guides .....	8
	Tool Options Guide.....	8
	Identifying Key Facial Features	
	The Landmark/Rotation Pictorial Guide.....	9
	The FAS Face and Other Anomalies Guide .....	9
	Measuring the Internal Measure of Scale	
	The Single-Distance Tool .....	10
	Measuring PFLs and Inner Canthal Distance	
	The 3-Distances Tool.....	10

Measuring Lip Thinness	
The Circularity Tool .....	11
The Lip-Philtrum Guides .....	11
Measuring Philtrum Smoothness	
The Lip-Philtrum Guides .....	12
Guide for Scoring the FAS Facial Phenotype	
The Facial ABC-Score and 4-Digit Rank Table .....	13
Measuring the Rotation of the Facial Image	
The Rotation Pictorial Guides.....	14
<b>PART V ANALYZING AN IMAGE</b>	
Example of a Complete Analysis and Outcome Report .....	15
Opening an Image and Worksheet	
To Open a New Photo.....	16
To Open an Existing Subject File .....	16
Enhancing the Image	
Brightness .....	17
Contrast .....	17
Rotation.....	17
Zoom.....	17
Clear overlays .....	17
Analyzing an Image	
Recording Data on the Worksheet .....	18
Patient Identification .....	18
Photograph Identification.....	19
Customizing Lists .....	20
Real Length of Internal Measure of Scale .....	20
Anomalies .....	20
Syndromes.....	20
Source of Image .....	20
Assessor of Image .....	20
Race of Subject .....	20
Racial Norms for Palpebral Fissure Lengths .....	20
Eyes.....	21
Lip/Philtrum.....	22
Other Anomalies / Syndromes .....	23
Photo Quality .....	24
Outcome.....	27
Practice Image Set “John Doe” .....	28
Circularity Practice Guide.....	29
<b>PART VI SAVING / DELETING / EXPORTING THE DATA</b>	
The Facial Software Database .....	30
Saving Data to the Database .....	31
Deleting Data from the Database.....	31
Exporting Data from the Database.....	32

**PART VII VIEWING / PRINTING THE OUTCOME REPORT**

Viewing the Outcome Report .....	33
Printing the Outcome Report .....	34
Outcome Report for Practice Case.....	35
<b>Dictionary</b> .....	36
<b>References</b> .....	39

Watch for Updates and Frequently Asked Questions on the FAS DPN Website

<http://depts.washington.edu/fasdpn>

## ACKNOWLEDGMENTS

---

The development of this software was supported in part by the following agencies and contributors:

Royalty Research Fund  
University of Washington  
Seattle Washington

Washington Research Foundation  
Seattle Washington

John B. Chavez FAS Fund

I would like to thank James M. Kinzel who served as the sole programmer for this software.

I would like to thank Kristen Daniels, M.L.I.S., Joshua Hunter, B.S. and Cliff Astley, M.S. for their invaluable technical assistance throughout the development of this program.

And special thanks are extended to all families who have supported and assisted the Washington State FAS Diagnostic & Prevention Network in our ongoing development and enhancement of FAS diagnostic and screening tools.