

**Scanning Probe Microscope
(Veeco Dimension 3100)
Training Notebook
(Advanced)**

***Nanotech User Facility
Center for Nanotechnology
March 2006***







Advanced Features

1. Motorized Sample Stage
2. Force Measurement
3. Point-and-Shoot
4. Q-Control
5. Fast Scan
6. Thermal Tune
7. Liquid Cell



1. Motorized Sample Stage

1. Set **Origin, X- and Y- Axis**.
 - Find a mark on the surface by moving the stage laterally using the arrows. Move the mark to the center of the crosshair.
 - Click on **Set Reference...** under Tools > Stage. Click on **Origin** to set the current position as the origin.
 - Move the stage horizontally to another mark and click on **X Axis** to set the X axis.
 - Bring the stage back to the Origin.
 - Click on **Move to (X, Y)** under Tools > Stage.
 - Type in 0 in **X Position** and **Y Position**.
 - Stage moves back to the Origin.
 - Move the stage vertically to another mark and click on **Y Axis** to set the Y axis.



1. Motorized Sample Stage (Cont.)

2. Move to a specific point on the sample.
 - Click on **Move to (X, Y)** under Tools > Stage.
 - Type in estimated values in **X Position** and **Y Position**, then click on **Go To**.
 - Type in values in **X Position Change** and **Y Position Change** to move the point to the center of the crosshair.
 - Start a quick scan (128 x128 pixels).
 - Use **Offset** to move the feature to the center of the image and click on **Execute**.
 - Adjust the **X position** and **Y Position** values according the Offset values and set the Offset values back to zero.
 - Rescan the image, the feature should be in the center.





Scan	
Scan size	500 nm
Aspect ratio	1.00
X offset	0.000 nm
Y offset	0.000 nm
Scan angle	0.00 °

1. Motorized Sample Stage (Cont.)

3. Programmed Move.

- Click on **Programmed Move** under Tools > Stage.
- Type in a new file name or browse an existing file.
- Add Step.
 - Click on **Add Step**. Step number changes.
 - Type in values in **X Position** and **Y Position**, then click on **Move**. Stage moves to the new position.
- Click on Add Step and then **Move To Selected**, **Step X** and **Step Y** values change.
- Repeat these steps to add new steps.
- **Remove Step**.
 - Type in the number of the step to be removed.
 - Click on Remove Step.
- Click on **Run** to run programmed moves.

(To capture images, click **Scan Parm List** window to active the tool bar. Use **Capture Continuous...**)



2. Force Measurement

1. Add **Ramp...**, **Ramp Parameters List...**, and **PicoForce...** to the workspace from the drop-down menu of **RealTime**.

2. Click on **Ramp** icon. 

3. Set **Data Type** and **X Data Type** in three channels as:
Defl. Vs. Time (or Z)
Z Sens. Vs. Time
Defl. Vs. Z Sens.




2. Force Measurement (Cont.)

Measure force on a single point once, or continuous.


1. Click on **Single**  or **Continuous**  force measurement icon.
Or, **Single** or **Continuous** in the panel.

2. Adjust tip travel distance (**Ramp Size**) and tip-surface interaction intensity (**Ramp Start**).



Ramp	
Ramp output	
Ramp size	949.3 nm
Z scan start	-398.1 nm



3. Click on **Stop** or  to stop the measurement.
4. Click on **Capture** to save force curves.



2. Force Measurement (Cont.)

Measure force on different points automatically.

1. Click on **Auto** tab to set up parameters.




2. Click on **Auto Force Measurement** icon to start the measurement.



2. Force Measurement (Cont.)

Measure forces on a single point according to the script.

1. Set **Strip Chart** parameters.
2. Click on **Script** tab to write script.





Mode	Value
Trigger mode	Off
Data type	Deflection
Trig threshold	0 V
Trig direction	Positive V
End mode	Retracted
Auto start	Disable
Surface delay	0.00 s
Retracted delay	0.00 s
Auto offset	Off
Strip chart rate	100 Hz
Strip chart size	10.0 s

(See Dimension Hybrid XYZ Scanning Probe Microscope Head manual, pp. 31 – 77, for detailed force measurement and analysis.)


3. Point-and-Shoot

1. Click **Point and Shoot** under RealTime.
2. Click on **Image/Scan** and **Capture Image** to pre-scan and save an image.
3. Click on **Point**, or **Line**, or **Box** tab to draw points, a line, or a box in the image. Click on **Convert To Points**.
4. Click **Capture** by Ramp for performing Ramp and Capture. (A subdirectory is created and the "before" image and force curves are saved in it.)
5. In the normal scan mode, scan and capture marked surface image after Point-and-Shoot.




3. Point-and-Shoot (Cont.)

6. Compare the point-and-shoot force curves with “after” image.
 - Rename “before” image.
 - Copy “after” image into the subdirectory and rename it as “before”.
7. Analyze forces.
 - Double click on folder icon.
 - Mark any force file.



(Dimension Hybrid XYZ Scanning Probe Microscope Head manual, pp. 40 – 52.)


4. Q Control

1. Click on **Tune** icon. 
2. Click on **Auto Tune**.
 - Click on **Compute Q** to get a Q value before Q control.



3. Set Q Control parameters.
4. Click on **More** and then **Adjust Q**.


Q CONTROL ACTIVE: Q = 10042 at 283.0 kHz with rms = 0.74 (Q factor = 33.44)



5. Turn off Q control before Engage.
 - More > Turn Off QC

(Engage is very slow or false engagement if QC is on.)
6. Click on Engage icon.

4. Q Control (Cont.)

7. Click on **Tune** icon  once Engage finishes.
8. In the pop-out window type in **1 μm** to lift the cantilever 1 μm away from the surface.
9. **Tune** window pops out. Click **More** > **Restore QC**, then **Exit**.
10. Scan images under Q control.

5. Fast scan




See NanoScope IVa Controller Manual, Software Version 6.13,
Section 7.4, pp. 61-72, for detail.

6. Thermal Tune

- It only works for cantilever resonance frequency less than 32 kHz. Usually works for soft cantilevers or in-liquid.
- See Atomic Force Microscope Thermal Tune Adapter manual for detail.

7. Liquid Cell

1. Choose a probe for scanning in liquid. (Si single crystal, FESP)
2. Mount the probe onto the **in-liquid probe holder**.
 - Put the holder on the **liquid cell docking station**.
 - Push the holder down to lift the spring
 - Insert the probe underneath the spring, align the probe against the back end of the groove.
3. Preparation.
 - Follow all steps in B.
(Note: More clearance between the bottom of the head and the sample stage is needed for in-liquid probe holder.)
4. Add liquid.
 - Take off the head.
 - Add a drop of liquid between the cantilever and the glass window of the in-liquid probe holder.
 - Add a drop of liquid onto the sample.
 - Put the head back to the dovetail trail.
(Note: The liquid should form a meniscus around the in-liquid probe holder.)



7. Liquid Cell (cont.)

5. Realignment.


- Realign laser to the free end of the cantilever.
- Readjust reflected laser to the center of the photo detector.
- Relocate the tip.

6. Contact Mode.

- Set Vertical Deflection to -1 V.
- Set Setpoint to 0 V.
- Set Integral Gain and Proportional Gain to 2.0 and 3.0.

7. Liquid Cell (cont.)

7. Tapping Mode.


- Set Vertical Deflection to 0 V, Integral Gain and Proportional Gain to 0.5 and 0.7.
- Tune cantilever *manually*.
- Click on Tune icon. 
- In the **Feedback Controls** panel, set **Z MODULATION** to DISABLED.
- Set **Sweep Control** parameters.

Drive Frequency	Sweep Width	Drive Amplitude	Amplitude Set Point	Amplitude Limit
10 kHz	20 kHz	2000 mV	0	2.5 V

- Use **Zoom In** and **Offset** to center the peak in the sweep plot.



- Adjust Drive Amplitude to get appropriate RMS Amplitude.
 - 0.3 -0.6 V for proteins.
 - 2 V for cells.



(Note: Two frequency ranges commonly used are 16-19 kHz and 8-12 kHz.)

7. Liquid Cell (cont.)

8. Engagement is the same as in-air.
9. Scan.
 - Adjust Setpoint to get better image.
 - Set Scan Rate slower in liquid.
10. Follow all steps in G. Shut Down.
11. Clean and dry the in-liquid probe holder.