

Zebrafish Neuromast Labeling Protocols

Nuclear Labels

| Dye | Ex_{max} | Em_{max} | Concentration | LM Filter | Laser Line | MRC-1024 Filter |
|------------|-------------------------|-------------------------|-----------------------|------------------|-------------------|------------------------|
| Syto 24/21 | 490/494 | 515/517 | 2 μ M – 6 μ M | FITC | 488 nm | 522DF35 |
| Yo-Pro-1 | 491 | 509 | 2 μ M – 6 μ M | FITC | 488 nm | 522DF35 |
| Lo-Pro-1 | 567 | 580 | 2 μ M – 6 μ M | Cy3 | 568 nm | 605DF32 |
| To-Pro-3 | 642 | 661 | 2 μ M – 6 μ M | Cy5 | 647 nm | 680DF32 |

Dyes are diluted in Embryo Medium from 1 mM stocks in DMSO. The stocks are kept frozen at -20° C. They may be aliquotted into small vials to avoid repeated freeze-thaw cycles and the attendant risks of water absorption and degradation. Neuromast sensory cells will stain within 20 minutes, although 1 hour seems better for timelapse imaging. The fluorescence of these dyes increases upon binding to DNA such that it is possible to use them without washing. All labels persist for at least 48 hours following incubation in live cells. Each of these nuclear labels are retained after fixation by 2% and 4% paraformaldehyde.

The monomeric cyanine dyes are cell impermeant in most cell types, which explains their apparent specificity for sensory cells. Their entry into live neuromast sensory cells is mediated by channels of unknown identity. The monomeric cyanine dyes display no increase in the number of cell types labeled by up to 2 hours of dye exposure. They can be employed as DNA counterstains in fixed samples. To-Pro-3 is particularly dependent upon concentrations over 2 μ M to resist bleaching. Lo-Pro is of limited utility due to a relatively broad spectra that fluoresces into the green channel upon excitation by blue wavelengths.

Syto dyes are permeant DNA labels that label nuclei only in live cells. Hair cells are strongly labeled and support cells are weakly labeled with Syto 24 at 30 min. Prolonged exposure to Syto 24, up to 2 hours, increases support cell nuclear intensity and labels additional cells in the epithelium. Syto dyes will not label DNA after aldehyde fixation. Syto dyes 21 and 24 were chosen over the other green Syto variants for desirable spectral properties and that they do not exhibit detectable RNA labeling.

Syto dyes are believed to interfere with transcription. Their effect on viability and development in long term live cell studies must be assessed for any experiment. Molecular Probes recommends that Syto dyes are only used in plastic containers to avoid non-specific binding to glass, and that phosphate-containing buffers be avoided for preparing working solutions.

Zebrafish Neuromast Labeling Protocols

Mitochondrial Labels

| Dye | Ex _{max} | Em _{max} | Concentration | LM Filter | Laser Line | MRC-1024 Filter |
|----------------------|-------------------|-------------------|---------------|--------------------|------------|---------------------|
| DASPEI | 461 | 589 | 5 drops/2 ml | Green LP or FM1-43 | 488 nm | 522DF35+ 605DF32 |
| Mitofluor Green | 490 | 516 | 25 nM | FITC | 488 nm | 522DF35 |
| Mitotracker Green FM | 490 | 516 | 25 nM | FITC | 488 nm | 522DF35 |
| Mitotracker Orange | 554 | 576 | 50 nM | Cy3 | 568 nm | 605DF32 |
| Mitotracker Red | 579 | 599 | 50 nM | Cy3/TR | 568 nm | 605DF32 |

DASPEI is prepared as a stock solution, 40 mg/100 ml distilled water, stored at 4° C. The Mitotracker and Mitofluor dyes are received as 50 µg vials to which anhydrous DMSO is added to produce a 1 mM stock. These stocks are kept frozen at -20°C. Small aliquots are recommended to avoid repeated freeze-thaw cycles. The reduced forms are particularly prone to oxidation. Working solutions of these dyes are dissolved in Embryo Medium. Label will appear in 15 minutes incubation, but they are commonly incubated for 1 hour when preparing fish for timelapse imaging.

DASPEI primarily labels hair cell mitochondria but also labels the nerve to the neuromast. The DASPEI emission curve is very broad, appearing in both “green” and “red” channels, giving yellow mitochondria and weakly red nuclei. Unlike the Mitotrackers, DASPEI will not label support cells and only labels a sub-population of hair cells in the neuromast. This label will leave the mitochondria after intense imaging on the confocal, presumably due to some form of phototoxic response causing mitochondrial depolarization.

The red and orange Mitotrackers will only label mitochondria possessing an electropotential difference relative to the cytoplasm. Mitotracker Green FM appears to accumulate regardless of mitochondrial membrane potential, making it a good marker for total mitochondrial mass. All Mitotrackers possess a weakly reactive chloromethyl group that forms thiol-peptide bonds with mitochondrial proteins. This allows them to be retained after fixation, however, they are also retained if the mitochondria become depolarized.

Mitofluors label mitochondria in a manner similar to the corresponding Mitotrackers. However, they are not retained by paraformaldehyde fixation. The Mitofluors are also retained upon depolarization of the mitochondria.

Higher concentrations of Mitotrackers and Mitofluor lead to dramatically increased non-specific background staining in surrounding epithelium as well as stain mitochondria in support cells and epithelial cells. The epithelial cells display a distinct mitochondrial

Zebrafish Neuromast Labeling Protocols

pattern that allows ready distinction from neuromast cells. The reduced versions, Red CM-H₂-XRos and Orange CM-H₂-TMRos, are preferred since they are non-fluorescent until entering an actively respiring mitochondrion and oxidized. Mitotracker Green FM is non-fluorescent until entering the lipid environment of the mitochondrial wall.

Transduction Channel Label

| Dye | Ex_{max} | Em_{max} | Concentration | LM Filter | Laser Line | MRC-1024 Filter |
|------------|-------------------------|-------------------------|----------------------|------------------|-------------------|------------------------|
| FM1-43 | 510 | 626 | 3 μM | FM1-43* | 488 nm | 605DF32 |
| FM1-43FX | 510 | 626 | 3 μM | Cy3 | 488 nm | 605DF32 |

*A Cy3 or Rhodamine filter will work. The Marianas FM1-43 filter uses a 490/20 excitation filter and the 617/73 emission filter. The Axioplan FM1-43 filter cube contains a 480/30 exciter and a 570 long pass filter.

FM1-43 is received as 100 μg dry powder reconstituted to 817 μM with 200 μl DMSO. It is dissolved in Embryo Medium for labeling zebrafish, with an incubation time of 30 to 45 seconds. Longer exposure to FM dyes may be toxic. While FM1-43 is not retained by paraformaldehyde fixation, FM1-43FX is a fixable version that displays identical excitation and emission spectra.

FM dyes are lipophilic labels generally absorbed by endocytosis with prolonged exposure. The internalized label then becomes generalized throughout the cell by membrane trafficking. The short incubation time used here insures that label is restricted to the neuromast sensory cells through a rapid uptake via the mechano-transduction channels. Membrane redistribution then transfers the dye throughout the cell over several minutes. The emission spectra for all FM dyes are extremely broad, a predominately red emission with a strong green signal, as well. It should also be noted that FM1-43 labeling seems to be inhibited by prior exposure to Syto dyes. This is being investigated further.

| Dye | Paraformaldehyde Fixable | Specificity |
|--------------------|---------------------------------|---------------------|
| Syto 24 | Yes | Support and Sensory |
| Yo-Pro-1 | Yes | Sensory |
| Lo-Pro-1 | Yes | Sensory |
| To-Pro-3 | Yes | Sensory |
| DASPEI | No | Mostly Sensory |
| Mitofluor Green | No | Sensory |
| Mitotracker Green | Yes | Sensory |
| Mitotracker Orange | Yes | Sensory |
| Mitotracker Red | Yes | Sensory |
| FM1-43 | No | Sensory |
| FM1-43FX | Yes | Sensory |

Zebrafish Neuromast Labeling Protocols

Working with fixed zebrafish.

Paraformaldehyde fixation of zebrafish induces strong autofluorescence in the green, red and far red channels. This autofluorescence may obscure all but exceptionally bright specific labels when viewed by widefield epifluorescence. Reducing the excitation intensity may assist in discernment of labels against background in fixed zebrafish when viewed with widefield epifluorescence. This autofluorescence can be reduced by using 2% paraformaldehyde or shorter fixation times.

High contrast, high resolution images from fixed fish require either confocal microscopy or deconvolution of image stacks collected from widefield epifluorescence.

The red Mitotrackers are based upon rhodamine derivatives. Anti-rhodamine antibodies from Molecular Probes are capable of labeling Mitotracker Red in zebrafish fixed with paraformaldehyde.