

**Stamatatos Lab**  
**Bulk Viral Infection in PBMCs**  
Contact: Zane Kraft at zane.kraft@seattlebiomed.org

**Day 1**

Isolate PBMCs and activate them with PHA for three days

**Day 3**

Count PBMCs

Spin the PBMCs at 2000rpm for 5 minutes and aspirate the supernatant

Resuspend cell pellet in RPMI w/ 20U/mL IL-2 at  $10 \times 10^6$  cells /mL

Aliquot 30 mLs of cells into a T225cm<sup>2</sup> flask

Add polybrene (2 ug/ml) to each tube

Incubate for 30 minutes at 37°C

Spin down the cells, 2000 rpm for 5 minutes, and aspirate the supernatant

Resuspend the cell pellet in 10mL of virus stock

Incubate at 37°C for 3 hours (gently mix cells every hour)

After 3 hours, spin down the cells at 2000 rpm for 5 minutes and aspirate the supernatant

Resuspend the cells in 100 mL of complete RPMI w/ 20 U/mL of IL-2

(this returns the cells to  $3 \times 10^6$  cells/mL)

Transfer to a T225cm<sup>2</sup> flask and incubate at 37 °C

**Day 7, 10, 14, and 17**

Take out 30-50 mL of bulk infection

Spin at 2000rpm for 5 minutes

Aliquot the supernatant into cryotubes (1 mL each) store at -80°C

Resuspend the cells in tube in 30-50 mL of complete RPMI w/ 20 U/mL IL-2 and return them to the flask

Repeat on days 10, 14 and 17 (Collect all supernatant on day 17)

NOTE: The viral inoculum used should only have been passaged once in PBMCs – primary isolates cannot be repeatedly passaged