

Important checks before starting a case

Its very important to make sure all system checks are working

Configuration Settings

Hardware Devices

Device	Status
Monitor (DataCaptor - localhost:7272)	Enabled

Status Watcher

Description	Status
Secure Server: uwmc-aims-ap.amc.uwm...	Connected
Monitor (DataCaptor - localhost:7272)	Connected
Message Handler	Ping Good

Service Watcher

Service	Location	Status
DSCache		Running
DSMessengerSvc		Running

Interaction Agents

Agent	Status
Decision Support Agent	Loaded
Primary Drug Event Detection Agent	Loaded
Patient-Specific ADE Detection Agent	Loaded
Tracked Item Detection Agent	Loaded
Bad Measurement Detection Agent	Loaded

[illegible]

Please verify all Status indicators are green before starting a case to insure proper communication between AIMS applications servers.

Any other color RED or YELLOW
Please notify Aims Support before
starting a case.

PH: 206.598.8695
PGR: 206.540.0937

All Items in the above image need to be green before starting a case. This is on the main merge startup page