

Installation Directions for the XP/VISTA compatible version of the FAS Facial Photographic Analysis Software

1. Installing the software:

- Login to XP/VISTA with an ADMINISTRATOR account
- Insert the CD with the Facial Software files if installing from a CD.
- Autorun will begin the process.
- If Autorun does not initiate the process, double click the Setup.exe file on the CD.
- If installing from the downloaded zip file instead, unzip and double click the Setup.exe file.
- Follow the onscreen instructions.
- You will likely encounter *Error 1911, Could not register type library for file C:\Program Files\Common files\Microsoft Shared\DAO\DAO350.DLL*. Click the Ignore button to continue with installation, and the software will install and run with no further problems.
- Click Finish.

The software is now loaded on your hard drive and a Facial Photographic Analysis folder will now appear on your desktop.

YOU MUST complete the next section before a non-administrator User can use the software. VISTA users must also complete the next section to use the software as a User or as an Administrator.

2. Altering User Permissions to enable the FAS software to write to its data folder in the system Program Files Folder.

Explanation. The FAS software saves the photo analysis data in an ACCESS file called data.fas that is located on C:\Program Files\FAS\Data. When the User uses the facial software, the software generates data that will be saved to data.fas. The software will also periodically compact this data.fas file. XP and VISTA do not permit programs run by a non-administrative User to “write” to folders in C:\Program Files, which means the facial software cannot modify its data.fas file. In addition, in VISTA, the facial software still is not able to modify the data.fas file in the FAS folder, even when run by an Administrator. To resolve these permission issues, individuals on both XP and VISTA systems must follow the instructions below.

Instructions for XP and VISTA. (Reminder: You must be an Administrator to do this).

- Double left click on the Computer icon.
- Double left click on the C drive.
- Double left click on the Program Files.
- Right click on the FAS folder. A menu will drop down.
- Left click on Properties. An “FAS Properties” window will open.
- Left click on the Security tab.
- In the Group and User Names: section, scroll down and left click Users.
- (VISTA users: Click the Edit button, then the Continue button to open the “Permissions for FAS” window and left click Users again.)
- Click the “Modify” checkbox in the Allow column.
- Click Apply
- Click OK (VISTA users: click OK again) to close the windows
- Close the computer window (click on the red x in the right upper corner).

To Open the Facial Software.

IMPORTANT! MS Access from the MS Office product must first be installed.

- Double click on the FAS Facial Analysis folder on your desktop.
- Double click on the FAS Facial Photographic Analysis Software file.
- A picture of the software will appear on your desktop.
- Wait a few seconds and the software will open to the main screen. Or you can click on the picture of the software to open it immediately.

Software Instruction Manual.

A copy of the Software Instruction Manual is in the FAS Facial Analysis Folder on your desktop.

Important Notice

READ BEFORE RE-INSTALLING THIS SOFTWARE – Reinstalling the software will overwrite your previous photo measurements saved in the data.fas database file created by the software. If you want to keep a copy of your previous photo measurements, you must follow the steps below **BEFORE REINSTALLING THE SOFTWARE**. If you not move and rename the data.fas file, the new software will overwrite that database file with a new (empty) database file. Follow these steps to move and rename the data.fas file.

- The database is called data.fas. It should be in the path C:\FAS\Data. If you do not see it there, simply search your hard drive for data.fas.
- Once you find it, rename the file (e.g., dataApril08.fas) and then move it to a location outside of C:\Program Files\FAS (e.g., place it on your desktop).
- Reinstall the software. When you do so, a new (empty) data.fas database will be created without overwriting the data you generated from previous photo analyses.
- Next, delete the newly created empty data.fas file created in C:\Program Files\FAS\Data. Move your renamed file (containing your data) to this location and rename it data.fas.

Frequently Asked Questions:

1. Is the software VISTA compatible?
Yes.
2. Can the software be used in a server environment?
No. The software was designed to be used on the computer it was loaded on.
3. Does the software work on a MacIntosh Computer?
No. The software requires Microsoft ACCESS, which is not supported on a MacIntosh Computer.