

**FAS Facial Photographic Analysis Software  
Version 1.0.0 (2003)**

**Reinstalling Instructions**

If you need to uninstall and reinstall the software and you want to save the data you have recorded, please follow the instructions below.

***BETA USERS NOTE:** The beta version of the database file will not work correctly with the upgraded software. Please start from scratch by uninstalling the beta version, deleting the Fetal Alcohol Syndrome file in C:\Program Files and installing the new version using the new data.fas file.*

- 1. Find the data.fas (database) file and move or copy it to somewhere else on your hard drive.**
  - a. Double click on the My Computer icon on your desktop.
  - b. Double click on LOCAL DISK (C:).
  - c. Double click on the Program Files folder.
  - d. Double click on the FAS folder.
  - e. Double click on the Data folder.
  - f. Highlight the data.fas file and move or copy it to somewhere else on your hard drive. It doesn't matter where, as long as it's not inside the FAS folder or any of the folders contained within that folder.
  - g. Make a note of where you put the data.fas file.
  
- 2. Uninstall the software.**
  - a. Click on the Start button on your toolbar.
  - b. Choose Settings.
  - c. Choose Control Panel
  - d. Double click on Add/Remove Programs.
  - e. Scroll down until you find FAS Facial Photographic Analysis Software.
  - f. Highlight the program name and click the Change/Remove button.
  - g. Follow the instructions in the uninstall wizard.
  
- 3. Install (or reinstall) the newer version of the software.**
  - a. Load the CD into your CD drive, and follow the instructions in your manual.
  
- 4. Replace the new data.fas file with the one you saved from the old installation.**
  - a. Go to wherever you put the old data.fas file. If you neglected to make a note of it, use Search in your Start menu to find it but remember that there will be two data.fas files now. You will want to find the one that is NOT located in C:\Program Files\FAS\Data.
  - b. Highlight the old data.fas file and move it to C:\Program Files\FAS\Data.
  - c. When asked, indicate that you'd like to replace the existing data.fas file.
  - d. Or if you prefer, you could delete the new data.fas file before moving the old one in.